



Product Information

LAUNDRY PRESOAK – Nappy Sanitiser

A premium quality oxygenated soak/wash sanitiser for pre wash treatment of stains and other contamination on linen and garments and for laundering and disinfecting soiled nappies. Laundry Presoak is **biodegradable** within internationally accepted criteria.

◆ DIRECTIONS:

For Linen:

Dissolve 25 grams of Laundry Presoak (approx slightly heaped tablespoon) in a $\frac{3}{4}$ full bucket of warm water. Add items & soak for 2-3 hours. Then remove and either hand or machine wash.

For Nappies:

Remove gross soil, then soak as above. Soak overnight if necessary. Rinse with water before machine washing.

Note:

Check garment label instructions before soaking. Not recommended for wool, silk or items with metal buttons or zippers.

◆ BIODEGRADABLE:

The surfactants in this product are biodegradable according to internationally accepted criteria.

◆ PRESENTATION:

Clark Laundry Presoak is available in:

- 2kg Tub
- 10kg pails
- 25kg bags

◆ HAZARDS:

Irritating to eyes. Repeated or prolonged contact may lead to skin irritation. Avoid inhaling product. Keep out of reach of children at all times.

◆ FIRST AID:

For advice, contact National Poisons Information Centre Phone 0800 764 766.

For eye contact: rinse cautiously with water for several minutes. Call a doctor if irritation persists.

For skin contact: wash with plenty of soap and water. If skin irritation persists get medical advice/attention.

A Safety Data Sheet is available for this product on request.